


Anti-Human Kappa Light Chain Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P01834

Item Number: CM0006768



Introduction

Rabbit monoclonal antibody against human kappa light chain constant region, validated for WB, IHC-P, IF-P, and ELISA. Ideal for immunoglobulin research, B cell studies, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Human Kappa Light Chain Rabbit mAb
Remarks/Alias	Km; HCAK1; IGKCD; Human Kappa Light Chain
Species	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human Kappa Light Chain (P01834)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01834
Protein Weight	12kDa
Shipping	Ice bag