

Anti-Human Lambda Light Chain Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P0CG04

Item Number: CM0006276



Introduction

Rabbit monoclonal antibody targeting human immunoglobulin lambda constant 1. Validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Ideal for adaptive immunity, B-cell, and light chain research.

[Learn More](#)

Parameter	Specification
Product Name	Human Lambda Light Chain Rabbit mAb
Species	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human Lambda Light Chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P0CG04
Protein Weight	11kDa
Shipping	Ice bag