

Anti-Human Lamp2 Monoclonal Antibody For Intracellular Flow Cytometry - P13473

Item Number: CM0008422



Introduction

This rabbit monoclonal anti-human LAMP2 antibody (conjugate) is validated for intracellular flow cytometry (FC). Targeting the lysosomal glycoprotein LAMP2 (CD107b), it is ideal for autophagy, lysosomal research, and immune presentation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LAMP2 mAb
Remarks/Alias	DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110
Species	Human
GeneID	3920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-375 of human LAMP2 (NP_002285.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P13473
Protein Weight	44kDa/45kDa
Shipping	Ice bag