

Anti-Human Lifr/Cd118 Monoclonal Antibody For Flow Cytometry - P42702

Item Number: CM0012100



Introduction

-conjugated rabbit monoclonal antibody targeting human LIFR/CD118, validated for flow cytometry. LIFR is a signal-transducing receptor for cytokines LIF, OSM, CNTF, CTF1, and CLCF1, involved in JAK/STAT pathway.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LIFR/CD118 mAb
Remarks/Alias	SWS; SJS2; STWS; CD118; LIF-R
Species	Human
GeneID (Human)	3977
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P42702
Protein Weight	124kDa
Shipping	Ice bag