

Anti-Human LILRA2/Cd85H Monoclonal Antibody For If/Icc, Fc - Q8N149

Item Number: CM0007984



Introduction

-conjugated rabbit monoclonal antibody targeting human LILRA2 (CD85h/ILT1/LIR7). Validated for immunofluorescence and flow cytometry. Ideal for innate immune response studies involving neutrophils and monocytes.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LILRA2/CD85h mAb
Remarks/Alias	ILT1; LIR7; CD85H; LIR-7
Species	Human
GeneID (Human)	11027
GeneID	11027
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-449 of human LILRA2/CD85h (NP_001124389.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q8N149
Protein Weight	53kDa/51kDa/48kDa/47kDa
Shipping	Ice bag