

Anti-Human Lilrb4/Cd85K/Ilt3 Monoclonal Antibody For Flow Cytometry - Q8Nhj6

Item Number: CM0010663



Introduction

Rabbit monoclonal antibody against human LILRB4/CD85k/ILT3 for flow cytometry applications. Detects the inhibitory immune receptor on monocytes, dendritic cells, and B cells. Suitable for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LILRB4/CD85k/ILT3 mAb
Remarks/Alias	B4; ILT3; LIR5; CD85K; ILT-3; LIR-5
Species	Human
GeneID (Human)	11006
GeneID	11006
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-259 of human LILRB4/CD85k/ILT3 (NP_001265355.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8NHJ6
Protein Weight	49kDa
Shipping	Ice bag