


Anti-Human Mmp9 Monoclonal Antibody For Intracellular Flow Cytometry - P14780

Item Number: CM0012420



Introduction

-conjugated rabbit anti-human MMP9 monoclonal antibody validated for intracellular flow cytometry. Targets matrix metalloproteinase-9 (gelatinase B), suitable for studies in extracellular matrix remodeling and leukocyte migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human MMP9 mAb
Remarks/Alias	MMP9; CLG4B; GELB; MANDP2; MMP-9; matrix metalloproteinase 9
Species	Human
Gene ID	4318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-707 of human MMP9 (NP_004985.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P14780
Protein Weight	78kDa
Shipping	Ice bag