


Anti-Human Mmp9 Rabbit Monoclonal Antibody For Flow Cytometry (Intracellular) - P14780

Item Number: CM0011620



Introduction

Recombinant rabbit monoclonal antibody targeting human matrix metalloproteinase-9 (MMP9). Validated for intracellular flow cytometry, it detects MMP9 involved in extracellular matrix degradation and leukocyte migration. Reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human MMP9 Monoclonal Antibody
Remarks/Alias	MMP9; CLG4B; GELB; MANDP2; MMP-9; matrix metalloproteinase 9
Species	Human
GeneID (Human)	4318
GeneID	4318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-707 of human MMP9 (NP_004985.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P14780
Protein Weight	78kDa
Shipping	Ice bag