

Anti-Human Nectin-4 Monoclonal Antibody () For Flow Cytometry - Q96Ny8

Item Number: CM0010368



Introduction

Rabbit monoclonal antibody against human Nectin-4, conjugated to . Validated for flow cytometry (FC). Suitable for studying cell adhesion, measles virus entry, and breast carcinoma research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Nectin-4 mAb
Remarks/Alias	LNIR; PRR4; EDSS1; PVRL4; nectin-4
Species	Human
Gene ID (Human)	81607
Gene ID	81607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-146 of human Nectin-4 (NP_112178.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96NY8
Protein Weight	24kDa/55kDa
Shipping	Ice bag