

Anti-Human Neurogranin (Nrgn) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92686

Item Number: CM0007498

Introduction

Rabbit monoclonal antibody targeting human Neurogranin (NRGN) for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for synaptic development and Alzheimer's disease research.

[Learn More](#)

Parameter	Value
Product Name	Neurogranin (NRGN) Rabbit mAb
Remarks/Alias	RC3; hng; Neurogranin (NRGN)
Species	Human
GeneID (Human)	4900
GeneID	4900
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-78 of human Neurogranin (NRGN) (NP_001119653.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92686
Protein Weight	8kDa
Shipping	Ice bag