

Anti-Human Neurogranin Rabbit Polyclonal Antibody For Wb, Elisa - Q92686

Item Number: CM0024595

Introduction

Rabbit polyclonal antibody against human neurogranin (NRGN/RC3) for WB and ELISA. Cross-reacts with mouse and rat. Calmodulin-binding protein involved in synaptic development and remodeling. Suitable for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Human Neurogranin Rabbit pAb
Remarks/Alias	RC3; hng; Human Neurogranin
Species	Human
Gene ID (Human)	4900
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 2-78 of human Human Neurogranin (NP_001119653.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92686
Protein Weight	8kDa
Shipping	Ice bag