


Anti-Human Parkin Polyclonal Antibody For Wb, If, Ihc, Elisa - O60260

Item Number: CM0018324



Introduction

High-quality rabbit polyclonal antibody against human Parkin, suitable for WB, IF, IHC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying parkin's role in mitophagy, Parkinson's disease, and ubiquitination pathways.

[Learn More](#)

Parameter	Value
Product Name	Human Parkin Rabbit pAb
Remarks/Alias	PDJ; AR-JP; LPRS2; PARK2; Parkin
Species	Human
GeneID (Human)	5071
GeneID	5071
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 81-280 of human Parkin (NP_004553.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60260
Protein Weight	52kDa
Shipping	Ice bag