


Anti-Human Parkin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60260

Item Number: CM0013191



Introduction

Rabbit polyclonal antibody targeting human Parkin (PRKN/PARK2), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for Parkinson's disease, mitophagy, and mitochondrial quality control research.

[Learn More](#)

Parameter	Value
Product Name	Human Parkin Rabbit pAb
Remarks/Alias	PDJ; AR-JP; LPRS2; PARK2; Parkin
Species	Human
GeneID (Human)	5071
GeneID	5071
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-387 of human Parkin (NP_004553.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60260
Protein Weight	52kDa
Shipping	Ice bag