

Anti-Human Pd-L1/Cd274 Monoclonal Antibody For Flow Cytometry - Q9Nzq7

Item Number: CM0009024



Introduction

Rabbit monoclonal antibody against human PD-L1/CD274, suitable for flow cytometry. Recognizes the 33 kDa programmed death-ligand 1, a key immune checkpoint protein (UniProt Q9NZQ7). Ideal for immuno-oncology research.

[Learn More](#)

Parameter	Details
Product Name	Anti-Human PD-L1/CD274 Monoclonal Antibody
Remarks/Alias	B7-H; B7H1; PDL1; PD-L1; hPD-L1; PDCD1L1; PDCD1LG1
Species	Human
GeneID (Human)	29126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PD-L1/CD274 (Q9NZQ7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NZQ7
Protein Weight	33kDa
Shipping	Ice bag