


Anti-Human Pla2G7/Paf-Ah/Lp-Pla2 Rabbit Polyclonal Antibody For Wb And Elisa - Q13093

Item Number: CM0023866



Introduction

Human PLA2G7/PAF-AH/Lp-PLA2 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse. Useful for studying inflammatory phospholipid metabolism and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	Human PLA2G7/PAF-AH/Lp-PLA2 Rabbit pAb
Remarks/Alias	PLA2G7; LDL-PLA2; LP-PLA2; PAFAD; PAFAH; phospholipase A2 group VII; Human PLA2G7/PAF-AH/Lp-PLA2
Species	Human
GeneID (human)	7941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-441 of human PLA2G7/PAF-AH/Lp-PLA2 (NP_005075.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13093
Protein Weight	50kDa
Shipping	Ice bag