

Anti-Human Podoplanin Monoclonal Antibody () For Flow Cytometry - Q86YI7

Item Number: CM0009228



Introduction

Rabbit monoclonal antibody against human podoplanin (PDPN), conjugated to , for flow cytometry. Reacts with human podoplanin; targets extracellular domain residues 21-123. Suitable for cell migration, adhesion, and EMT research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Podoplanin mAb
Remarks/Alias	PDPN; AGGRUS; GP36; GP40; Gp38; HT1A-1; OTS8; PA2.26; T1A; T1A-2; T1A2; T11A; podoplanin
Species	Human
GeneID (human)	10630
GeneID	10630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-123 of Human Podoplanin (NM_006474.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q86YL7
Protein Weight	12kDa/16kDa/18kDa/24kDa
Shipping	Ice bag