

Anti-Human Pro-Collagen I Alpha (Col1A1) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P - P02452

Item Number: CM0020355



Introduction

Reliable rabbit polyclonal antibody against human pro-collagen I alpha (COL1A1), suitable for WB, IF-P, and IHC-P. Designed for studying collagen type I in extracellular matrix, bone, and connective tissue biology.

[Learn More](#)

Parameter	Value
Product Name	Human Pro-Collagen I alpha Rabbit pAb
Remarks/Alias	COL1A1; EDSC; OI1; OI2; OI3; OI4; collagen alpha-1(I) chain; Human Pro-Collagen I alpha
Species	Human
Gene ID	1277
Immunogen	Recombinant fusion protein corresponding to amino acids 26-161 of human Collagen I/COL1A1 (NP_000079.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P
Cross Reactivity	Human
SWISS	P02452
Protein Weight	138kDa
Shipping	Ice bag