


Abflo 488 Anti-Human Rankl/Cd254 Monoclonal Antibody For Flow Cytometry (Intracellular) - O14788

Item Number: CM0009276



Introduction

ABflo 488-conjugated rabbit monoclonal antibody targeting human RANKL/CD254 (UniProt: O14788). Validated for intracellular flow cytometry. Ideal for investigating osteoclastogenesis, bone resorption, and immune cell crosstalk.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human RANKL/CD254 mAb
Remarks/Alias	ODF; OPGL; sOdf; CD254; OPTB2; RANKL; TNLG6B; TRANCE; hRANKL2
Species	Human
GeneID	8600
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 136-317 of human RANKL/CD254 (NP_003692.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (Intracellular)
Cross Reactivity	Human
SWISS	O14788
Protein Weight	27kDa/30kDa/35kDa
Shipping	Ice bag