

Anti-Human Siglec9 Monoclonal Antibody For Flow Cytometry - Q9Y336

Item Number: CM0010993

Introduction

-conjugated rabbit anti-human SIGLEC9 (CD329) monoclonal antibody optimized for flow cytometry. Detects human SIGLEC9, a sialic acid-binding lectin with roles in cell adhesion and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human SIGLEC9 mAb
Remarks/Alias	CD329; CDw329; FOAP-9; siglec-9; OBBP-LIKE
Species	Human
GenelD	27180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-348 of human SIGLEC9 (NP_055256.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9Y336
Protein Weight	50kDa, 52kDa
Shipping	Ice bag