

Anti-Human Syt12 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8lv01

Item Number: CM0019390

Introduction

Anti-Human SYT12 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects human SYT12 (47 kDa), cross-reacts with mouse. Ideal for synaptic transmission studies.

[Learn More](#)

Parameter	Value
Product Name	SYT12 Rabbit pAb
Remarks/Alias	SYT11; sytXII; SYT12
Species	Human
GeneID	91683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 95-200 of human SYT12 (NP_001171351.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q8IV01
Protein Weight	47kDa
Shipping	Ice bag