


Anti-Human Tim-3/Havcr2 Monoclonal Antibody For Fc - Q8Tdq0

Item Number: CM0007853



Introduction

Rabbit monoclonal antibody against human TIM-3/HAVCR2, validated for flow cytometry.

Suitable for immunology research, targeting T-cell immunoglobulin mucin receptor 3 (CD366) involved in immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TIM-3/HAVCR2 mAb
Remarks/Alias	TIM3; CD366; KIM-3; SPTCL; TIMD3; Tim-3; TIMD-3; HAVcr-2
Species	Human
GeneID (Human)	84868
GeneID	84868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-200 of human TIM-3/HAVCR2 (NP_116171.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TDQ0
Protein Weight	33kDa
Shipping	Ice bag