


Anti-Human Tnfrsf11A Polyclonal Antibody For Wb, Elisa - Q9Y6Q6

Item Number: CM0023000



Introduction

High-quality rabbit polyclonal antibody against human TNFRSF11A (RANK) for WB and ELISA. Cross-reacts with mouse and rat. Ideal for osteoclastogenesis and bone research applications.

[Learn More](#)

Parameter	Value
Product Name	Human TNFRSF11A Rabbit pAb
Remarks/Alias	FEO; OFE; ODFR; OSTs; PDB2; RANK; CD265; OPTB7; TRANCER; LOH18CR1; TRANCE-R; Human TNFRSF11A
Species	Human
Gene ID (Human)	8792
Gene ID	8792
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-212 of human TNFRSF11A (NP_003830.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6Q6
Protein Weight	66kDa
Shipping	Ice bag