

Anti-Human Trem1 (Cd354) Monoclonal Antibody For Flow Cytometry - Q9Np99

Item Number: CM0009020



Introduction

-conjugated rabbit monoclonal antibody targeting human TREM1 (CD354), a myeloid receptor that amplifies inflammatory responses. Designed for flow cytometry. Suitable for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TREM1 mAb
Remarks/Alias	CD354; TREM-1
Species	Human
Gene ID (Human)	54210
Gene ID	54210
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-200 of human TREM1 (NP_061113.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NP99
Protein Weight	18kDa/26kDa
Shipping	Ice bag