


Anti-Human Trka Monoclonal Antibody () For Flow Cytometry - P04629

Item Number: CM0008589



Introduction

-conjugated rabbit anti-human TrkA monoclonal antibody validated for flow cytometry on human samples. Targets the high affinity nerve growth factor receptor, ideal for neurobiology and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human TrkA mAb
Remarks/Alias	NTRK1; MTC; TRK; TRK1; TRKA; Trk-A; p140-TrkA; neurotrophic receptor tyrosine kinase 1
Species	Human
GeneID	4914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-407 of human TrkA (NP_002520.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04629
Protein Weight	77kDa/83kDa/86kDa/87kDa
Shipping	Ice bag