

Anti-Human Vcam-1/Cd106 Monoclonal Antibody For Flow Cytometry - P19320

Item Number: CM0011908



Introduction

-conjugated rabbit monoclonal antibody targeting human VCAM-1/CD106, validated for flow cytometry (FC). Ideal for studying leukocyte adhesion, endothelial inflammation, and immune surveillance. Recognizes human VCAM-1 protein (P19320).

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human VCAM-1/CD106 mAb
Remarks/Alias	CD106
Species	Human
Gene ID (Human)	7412
Gene ID	7412
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-698 of human VCAM1 (NP_001069.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P19320
Protein Weight	71kDa/74kDa/81kDa
Shipping	Ice bag