

Anti-Human/Monkey Cd3 Mab For Flow Cytometry - P07766

Item Number: CM0008931



Introduction

-conjugated rabbit anti-human/monkey CD3 monoclonal antibody for flow cytometry detection of CD3 epsilon on T cells. Validated for human and cynomolgus monkey samples. Ideal for immunophenotyping studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD3 mAb
Remarks/Alias	T3E, TCRE, IMD18, CD3epsilon
Species	Human
GeneID (Human)	916
GeneID	916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human CD3 (NP_000724.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag