

Anti-Human/Monkey Cd64 Monoclonal Antibody For Flow Cytometry - P12314

Item Number: CM0010309



Introduction

-conjugated rabbit monoclonal antibody targeting human/monkey CD64 (FCGR1A) for flow cytometry. Recognizes the high-affinity IgG receptor expressed on monocytes/macrophages. SWISS: P12314.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD64 mAb
Remarks/Alias	FCGR1A; CD64; CD64A; FCRI; IGFR1; Fc fragment of IgG receptor Ia
Species	Human
GeneID (Human)	2209
GeneID	2209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-288 of human CD64 (NP_000557.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P12314
Protein Weight	38kDa/42kDa
Shipping	Ice bag