

Anti-Human/Monkey Cd69 Monoclonal Antibody For Flow Cytometry - Q07108

Item Number: CM0009814



Introduction

Rabbit monoclonal antibody against CD69, directly conjugated with . Specific for human and cynomolgus CD69. Validated for flow cytometry to study early activation of T cells, NK cells, and other leukocytes.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD69 mAb
Remarks/Alias	AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26
Species	Human
GeneID (Human)	969
GeneID	969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 78-199 of human CD69 (NP_001772.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q07108
Protein Weight	23kDa
Shipping	Ice bag