


Anti-Human/Monkey Cd8A Monoclonal Antibody For Flow Cytometry - P01732

Item Number: CM0011892



Introduction

Rabbit monoclonal antibody targeting human and cynomolgus monkey CD8a, validated for flow cytometry. Recognizes the CD8 alpha chain, a coreceptor for MHC class I, essential for T-cell-mediated immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD8a mAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha
Species	Human
Gene ID	925 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8A (NP_001759.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01732
Protein Weight	21kDa/25kDa/30kDa
Shipping	Ice bag