


Anti-Human/Monkey IgD Monoclonal Antibody () For Flow Cytometry - P01880

Item Number: CM0012564



Introduction

High-quality -conjugated rabbit anti-human/monkey IgD monoclonal antibody validated for flow cytometry (FC). Recognizes human and cynomolgus IgD, ~47 kDa. Ideal for B-cell immunophenotyping and humoral immunity studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgD mAb
Species	Human
Gene ID	3495
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01880
Protein Weight	47kDa
Shipping	Ice bag