

Anti-Human/Monkey IgD Monoclonal Antibody For Flow Cytometry - P01880

Item Number: CM0011036



Introduction

-conjugated rabbit monoclonal anti-human/monkey IgD antibody for flow cytometry. Detects the delta heavy chain constant region; useful for B-cell immunophenotyping and humoral immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgD mAb
Remarks/Alias	IGHD
Species	Human
GeneID	3495
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01880
Protein Weight	47kDa
Shipping	Ice bag