

Anti-Human/Mouse Cd34 Rabbit Polymab® Monoclonal Antibody For Wb, Ihc-P, Elisa - P28906 / Q64314

Item Number: CM0007706

Introduction

Rabbit PolymAb® monoclonal antibody targeting human/mouse CD34 antigen. Validated for WB, IHC-P, ELISA. Detects a 41 kDa band. Ideal for hematopoietic stem cell, adhesion, and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	CD34 Rabbit PolymAb®
Remarks/Alias	CD34;CD34 molecule;GIG3;MORT1
Species	Human/Mouse
Gene ID (Human)	947
Gene ID	947/12490
Immunogen	Recombinant protein of human CD34
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28906/Q64314
Protein Weight	41kDa
Shipping	Ice bag