

Anti-Huntingtin Monoclonal Antibody For Wb, If-F, Elisa - P42858

Item Number: CM0002973



Introduction

Rabbit monoclonal antibody targeting Huntingtin, validated for Western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat HTT. Ideal for studying Huntington's disease and autophagy pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Huntingtin Rabbit mAb
Remarks/Alias	HD; IT15; LOMARS; Huntingtin
Species	Human
GeneID	3064
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Huntingtin (P42858)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42858
Protein Weight	348 kDa
Shipping	Ice bag