

Anti-Huntingtin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42858

Item Number: CM0022755



Introduction

Rabbit polyclonal antibody against human Huntingtin (aa 435-635). Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. For research use in Huntington's disease and neurodegenerative studies.

[Learn More](#)

Parameter	Value
Product Name	Huntingtin Rabbit pAb
Remarks/Alias	HD; IT15; LOMARS; Huntingtin
Species	Human
Gene ID (Human)	3064
Gene ID	3064
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 435-635 of human Huntingtin (NP_002102.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42858
Protein Weight	348kDa
Shipping	Ice bag