

Anti-Hur/Elavl1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q15717

Item Number: CM0003077



Introduction

High-quality rabbit monoclonal antibody targeting HuR/ELAVL1, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Recognizes human, mouse, and rat ELAVL1. Suitable for mRNA stability and post-transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HuR/ELAVL1 Rabbit mAb
Remarks/Alias	HUR; Hua; MeIG; ELAV1; L1
Species	Human
GeneID (Human)	1994
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HuR/ELAVL1 (Q15717).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15717
Protein Weight	36kDa
Shipping	Ice bag