


Anti-Hur/Elavl1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15717

Item Number: CM0030548



Introduction

Rabbit polyclonal HuR/ELAVL1 antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human and mouse HuR, an RNA-binding protein that stabilizes AU-rich element-containing mRNAs for gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HuR/ELAVL1 Rabbit pAb
Remarks/Alias	HUR; Hua; MeIG; ELAV1; L1
Species	Human
Gene ID (Human)	1994
Gene ID	1994
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HuR/ELAVL1 (NP_001410.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15717
Protein Weight	36kDa
Shipping	Ice bag