


Anti-Hur/Elavl1 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Elisa - Q15717

Item Number: CM0024275



Introduction

KO validated HuR/ELAVL1 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Recognizes RNA-binding protein HuR, stabilizes ARE-containing mRNAs; ideal for post-transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HuR/ELAVL1 Rabbit pAb
Remarks/Alias	HUR; Hua; MeIG; ELAV1; L1
Species	Human
GeneID (Human)	1994
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human HuR/ELAVL1 (NP_001410.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15717
Protein Weight	36kDa
Shipping	Ice bag