


Anti-Hus1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O60921

Item Number: CM0028139



Introduction

Rabbit polyclonal antibody targeting human HUS1 (O60921), a 32kDa checkpoint protein. Validated for WB, IF/ICC, IP, ELISA; cross-reacts with human and mouse. Ideal for DNA repair and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	HUS1 Rabbit pAb
Remarks/Alias	hHUS1; HUS1
Species	Human
GeneID (human)	3364
GeneID	3364
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human HUS1 (NP_004498.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60921
Protein Weight	32kDa
Shipping	Ice bag