


Anti-Hus1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60921

Item Number: CM0014131



Introduction

Rabbit polyclonal antibody against human HUS1 checkpoint protein (32kDa). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for DNA damage and cell cycle checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	HUS1 Rabbit pAb
Remarks/Alias	hHUS1; HUS1
Species	Human
Gene ID	3364
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human HUS1 (NP_004498.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60921
Protein Weight	32kDa
Shipping	Ice bag