

Anti-Hus1B Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nhy5

Item Number: CM0016946

Introduction

Rabbit polyclonal antibody targeting human checkpoint protein HUS1B, validated for western blot (WB) and ELISA. Demonstrates cross-reactivity with mouse and rat, making it suitable for multi-species studies of cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	HUS1B Rabbit pAb
Remarks/Alias	HUS1B
Species	Human
GeneID (Human)	135458
GeneID	135458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-278 of human HUS1B (NP_683762.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8NHY5
Protein Weight	31kDa
Shipping	Ice bag