


Anti-Huwe1 Polyclonal Antibody For Wb And Elisa - Q7Z6Z7

Item Number: CM0019259



Introduction

Rabbit polyclonal antibody against human HUWE1, validated for Western blot and ELISA. Detects HUWE1 (~482kDa) in human and rat samples. Ideal for ubiquitination and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HUWE1 Rabbit pAb
Remarks/Alias	MULE; Ib772; LASU1; MRXST; UREB1; HECTH9; URE-B1; ARF-BP1; HSPC272; HUWE1
Species	Human
Gene ID	10075
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HUWE1 (NP_113584.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q7Z6Z7
Protein Weight	482kDa
Shipping	Ice bag