

Anti-Hvem/Tnfrsf14 Rabbit Pab For Wb, Elisa - Q92956

Item Number: CM0020975



Introduction

Rabbit polyclonal antibody against human HVEM/TNFRSF14, validated for Western blot and ELISA. This 30 kDa membrane receptor mediates immune costimulation, antiviral responses, and herpesvirus entry. Suitable for immunology, infectious disease, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	HVEM/TNFRSF14 Rabbit pAb
Remarks/Alias	TR2; ATAR; HVEA; HVEM; CD270; LIGHTR; HVEM/TNFRSF14
Species	Human
GeneID (Human)	8764
GeneID	8764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-202 of human HVEM/TNFRSF14 (NP_003811.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92956
Protein Weight	30kDa
Shipping	Ice bag