

Anti-Hyal1 Monoclonal Antibody For Wb, Elisa - Q12794

Item Number: CM0004369



Introduction

Rabbit monoclonal HYAL1 antibody for Western blot (WB) and ELISA, reactive with human. Targets hyaluronidase-1, involved in tumor progression and lysosomal function. Suitable for cancer and mucopolysaccharidosis research.

[Learn More](#)

Parameter	Value
Product Name	HYAL1 Rabbit mAb
Remarks/Alias	HYAL1; HYAL-1; LUCA1; MPS9; NAT6; hyaluronidase-1
Species	human
Gene ID (Human)	3373
Gene ID	3373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-435 of human HYAL1 (NP_149349.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q12794
Protein Weight	10kDa/19kDa/23kDa/27kDa/37kDa/45kDa/48kDa
Shipping	Ice bag