

## Anti-Iah1 Rabbit Polyclonal Antibody For Wb, Elisa - Q2Taa2

Item Number: CM0013947

### Introduction

Rabbit polyclonal antibody targeting IAH1, a probable lipase involved in lipid metabolism. Validated for WB and ELISA applications with cross-reactivity to human, mouse, and rat. Detects a 28kDa band. Ideal for hydrolase and lipid degradation research.

[Learn More](#)

| Parameter        | Value                 |
|------------------|-----------------------|
| Product Name     | IAH1 Rabbit pAb       |
| Remarks/Alias    | IAH1                  |
| Species          | Human                 |
| GeneID (Human)   | 285148                |
| GeneID           | 285148                |
| Immunogen        | Recombinant protein   |
| Source           | Rabbit                |
| Category         | Polyclonal Antibodies |
| Application      | WB, ELISA             |
| Cross Reactivity | Human, Mouse, Rat     |
| SWISS            | Q2TAA2                |
| Protein Weight   | 28kDa                 |
| Shipping         | Ice bag               |