

Anti-Iapp (Amylin, Dap) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P10997

Item Number: CM0027580



Introduction

Rabbit polyclonal antibody against human IAPP (Amylin/DAP), suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Useful for studies on energy homeostasis and type 2 diabetes amyloid aggregation.

[Learn More](#)

Parameter	Value
Product Name	IAPP Rabbit pAb
Remarks/Alias	DAP; IAP; IAPP
Species	Human
Gene ID (Human)	3375
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-89 of human IAPP (NP_000406.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10997
Protein Weight	10 kDa
Shipping	Ice bag