

Anti-Ibsp Rabbit Polyclonal Antibody For Wb, Elisa - P21815

Item Number: CM0020543

Introduction

Rabbit polyclonal antibody targeting human integrin-binding sialoprotein (IBSP), validated for Western blotting and ELISA, with cross-reactivity to mouse and rat. Suitable for bone mineralization and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	IBSP Rabbit pAb
Remarks/Alias	BSP; BNSP; SP-II; BSP-II; IBSP
Species	Human
GeneID (Human)	3381
GeneID	3381
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human IBSP (NP_004958.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21815
Protein Weight	35kDa
Shipping	Ice bag