

Anti-Ibtk Rabbit Polyclonal Antibody For Wb And Elisa - Q9P2D0

Item Number: CM0015886

Introduction

Anti-IBTK rabbit polyclonal antibody validated for WB and ELISA. Targets human IBTK (151kDa), an inhibitor of BTK kinase involved in B-cell development and NF-κB transcription.

[Learn More](#)

Parameter	Value
Product Name	IBTK Rabbit pAb
Remarks/Alias	BTKI; BTBD26; IBTK
Species	Human
GeneID (human)	25998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1200 of human IBTK (NP_056340.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9P2D0
Protein Weight	151kDa
Shipping	Ice bag