

## Anti-Icam-1/Cd54 Monoclonal Antibody For Flow Cytometry - Q00238

Item Number: CM0011329

### Introduction

Rabbit monoclonal antibody against rat ICAM-1/CD54, conjugated to and validated for flow cytometry. Ideal for studying leukocyte adhesion, inflammation, and immune cell trafficking.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Rat ICAM-1/CD54 mAb                                                                                       |
| Remarks/Alias    | CD54; ICAM; Icam1                                                                                                     |
| Species          | Rat                                                                                                                   |
| Gene ID (Human)  | 3383                                                                                                                  |
| Gene ID          | 25464                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 28-493 of rat ICAM-1/CD54 (NP_037099.1) |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | FC                                                                                                                    |
| Cross Reactivity | Rat                                                                                                                   |
| SWISS            | Q00238                                                                                                                |
| Protein Weight   | 60kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |