


Anti-Icam-1/Cd54 Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P05362 / P13597

Item Number: CM0011215



Introduction

Recombinant rabbit monoclonal antibody against ICAM-1 (CD54) for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Suited for leukocyte adhesion, inflammation, and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	ICAM-1/CD54 Rabbit PolymAb®
Remarks/Alias	MAL; CD54; Icam; Ly-4; Ly-47; Icam-1; MALA-2
Species	Human/Mouse
Gene ID (Human)	3383
Gene ID	3383/15894
Immunogen	Recombinant protein of human/mouse ICAM-1/CD54.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05362/P13597
Protein Weight	58kDa/59kDa
Shipping	Ice bag