


Anti-Icam-1/Cd54 Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - P05362

Item Number: CM0002432



Introduction

High-specificity recombinant rabbit monoclonal antibody targeting human ICAM-1/CD54 (UniProt P05362). Validated for WB, IHC-P, FC, and ELISA. Ideal for studying leukocyte adhesion, endothelial activation, and viral entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	ICAM-1/CD54 Rabbit mAb
Remarks/Alias	BB2; CD54; P3.58; MALA2; MyD10; ICAM-1/CD54
Species	Human
Gene ID (Human)	3383
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ICAM-1/CD54 (P05362))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P05362
Protein Weight	58kDa
Shipping	Ice bag